


Smac/Diablo Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Nr28

Item Number: CM0007147



Introduction

Smac/Diablo Rabbit mAb (Q9NR28) validated for WB, IF/ICC, and ELISA. Targets human, mouse, rat Diablo protein (27kDa). Ideal for apoptosis and IAP pathway research.

[Learn More](#)

Parameter	Value
Product Name	Smac/Diablo Rabbit mAb
Remarks/Alias	SMAC; DFNA64; Smac/Diablo
Species	Human
Gene ID (Human)	56616
Gene ID	56616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 56-239 of human Smac/Diablo (NP_063940.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR28
Protein Weight	27kDa
Shipping	Ice bag