


Smac/Diablo Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nr28

Item Number: CM0010538



Introduction

Recombinant rabbit monoclonal antibody against human Smac/Diablo (SMAC), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Smac/Diablo Rabbit mAb
Remarks/Alias	SMAC; DFNA64; Smac/Diablo
Species	Human
Gene ID (Human)	56616
Gene ID	56616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 56-239 of human Smac/Diablo (NP_063940.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR28
Protein Weight	27kDa/21kDa/22kDa
Shipping	Ice bag