


Smad2 Rabbit Pab - Q15796

Item Number: CM0017119



Introduction

Rabbit polyclonal antibody specific for human Smad2, suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse. Validated for research use in TGF-beta signaling and colorectal cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Smad2 Rabbit pAb
Remarks/Alias	JV18, LDS6, CHTD8, MADH2, MADR2, JV18-1, hMAD-2, hSMAD2, Smad2
Species	Human
GeneID (Human)	4087
GeneID	4087
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-270 of human Smad2 (NP_005892.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15796
Protein Weight	52kDa
Shipping	Ice bag