


Smad5 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q99717

Item Number: CM0003703



Introduction

Smad5 Rabbit Monoclonal Antibody for WB, IP, ELISA. Recognizes human and mouse Smad5. Smad5 is a transcriptional mediator in BMP signaling, critical for embryonic development and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	Smad5 Rabbit mAb
Remarks/Alias	DWFC; JV5-1; MADH5; Smad5
Species	Human
GeneID	4090
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human Smad5 (Q99717).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99717
Protein Weight	52kDa
Shipping	Ice bag