


Smc4 Rabbit Pab For Wb, Elisa - Q9Ntj3

Item Number: CM0017311



Introduction

SMC4 rabbit polyclonal antibody for WB and ELISA, targeting human SMC4 (CAP-C, SMC4L1) central to chromosome condensation. Uniprot Q9NTJ3, predicted MW 147kDa.

[Learn More](#)

Parameter	Value
Product Name	SMC4 Rabbit pAb
Remarks/Alias	CAPC; CAP-C; SMC-4; SMC4L1; SMC4
Species	Human
GeneID (Human)	10051
GeneID	10051
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SMC4 (NP_001275682.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NTJ3
Protein Weight	147kDa
Shipping	Ice bag