


Smg5 Rabbit Pab - Q9Upr3

Item Number: CM0021875



Introduction

Rabbit polyclonal antibody against human SMG5, validated for Western Blot and ELISA. Cross-reacts with mouse and rat. Suitable for nonsense-mediated mRNA decay research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	SMG5 Rabbit pAb
Remarks/Alias	EST1B; SMG-5; LPTSRP1; LPTS-RP1; SMG5
Species	Human
GeneID (Human)	23381
GeneID	23381
Immunogen	Synthetic peptide (a synthetic peptide corresponding to a sequence within amino acids 950-1016 of human SMG5 (NP_056142.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UPR3
Protein Weight	114kDa
Shipping	Ice bag