


Smn2 Rabbit Pab - Q16637

Item Number: CM0023116



Introduction

Rabbit polyclonal antibody to SMN2 (survival motor neuron protein), validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Ideal for studying snRNP assembly and splicing in neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	SMN2 Rabbit pAb
Remarks/Alias	SMNC; BCD541; GEMIN1; TDRD16B; C-BCD541; SMN2
Species	Human
Gene ID (Human)	6607
Gene ID	6607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human SMN2 (NP_059107.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16637
Protein Weight	32kDa
Shipping	Ice bag