

Smndc1 Rabbit Polyclonal Antibody For Wb, Elisa - O75940

Item Number: CM0014817



Introduction

SMNDC1 Rabbit Polyclonal Antibody validated for Western blotting and ELISA, detecting human, mouse, and rat SMNDC1 protein. Ideal for splicing research and biomarker studies.

[Learn More](#)

Parameter	Details
Product Name	SMNDC1 Rabbit pAb
Remarks/Alias	SMNR; SPF30; TDRD16C; SMNDC1
Species	Human
GeneID (Human)	10285
GeneID	10285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of human SMNDC1 (NP_005862.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75940
Protein Weight	27kDa
Shipping	Ice bag