

Smpd4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nxe4

Item Number: CM0013077

Introduction

SMPD4 Rabbit Polyclonal Antibody (pAb) targeting human SMPD4 (SKNY, nSMase-3). Validated for WB, IF/ICC, ELISA with cross-reactivity to human and mouse. Ideal for sphingolipid metabolism and membrane homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	SMPD4 Rabbit pAb
Remarks/Alias	SKNY; NET13; NEDMABA; NEDMEBA; NSMASE3; NSMASE-3; SMPD4
Species	Human
Gene ID (Human)	55627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-800 of human SMPD4 (NP_060421.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NXE4
Protein Weight	98kDa
Shipping	Ice bag