

Smyd1 Rabbit Pab - Q8Nb12

Item Number: CM0024944

Introduction

SMYD1 Rabbit pAb (WB, ELISA) detects human, mouse, and rat SMYD1. This polyclonal antibody targets the histone-lysine N-methyltransferase SMYD1, a transcriptional repressor essential for cardiomyocyte differentiation and cardiac morphogenesis.

[Learn More](#)

Parameter	Specification
Product Name	SMYD1 Rabbit pAb
Remarks/Alias	BOP; KMT3D; ZMYND18; ZMYND22; SMYD1
Species	Human
Gene ID (Human)	150572
Gene ID	150572
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NB12
Protein Weight	57kDa
Shipping	Ice bag