

Smyd4 Rabbit Pab - Q8Iyr2

Item Number: CM0022256

Introduction

Rabbit polyclonal antibody against human SMYD4 (ZMYND21). Suitable for WB, IF/ICC, IP, ELISA. Recognizes human and mouse SMYD4. Ideal for cardiac development and epigenetic regulation research.

[Learn More](#)

Parameter	Value
Product Name	SMYD4 Rabbit pAb
Remarks/Alias	ZMYND21; SMYD4
Species	Human
Gene ID (Human)	114826
Gene ID	114826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human SMYD4 (NP_443160.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IYR2
Protein Weight	89kDa
Shipping	Ice bag