

Smyd4 Rabbit Polyclonal Antibody For Wb - Q8Iyr2

Item Number: CM0029434

Introduction

Rabbit polyclonal antibody targeting SMYD4, a protein-lysine N-methyltransferase involved in cardiac development. Validated for WB, IHC-P, IF/ICC, IP, ELISA in Human, Mouse, Rat. Ideal for epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	SMYD4 Rabbit pAb
Remarks/Alias	ZMYND21; SMYD4
Species	Human
Gene ID (Human)	114826
Gene ID	114826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human SMYD4 (NP_443160.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IYR2
Protein Weight	89kDa
Shipping	Ice bag