

Snai3 Rabbit Polyclonal Antibody For Wb And Elisa - Q3Knw1

Item Number: CM0022183

Introduction

High-quality rabbit polyclonal antibody against human SNAI3, validated for Western blotting and ELISA. Suitable for human, mouse, and rat samples. Ideal for studying transcriptional repression and myoblast differentiation.

[Learn More](#)

Parameter	Value
Product Name	SNAI3 Rabbit pAb
Remarks/Alias	SMUC; SNAI3; ZNF293; Zfp293; SNAI3
Species	Human
GeneID (Human)	333929
GeneID	333929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 7-186 of human SNAI3 (NP_840101.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q3KNW1
Protein Weight	32kDa
Shipping	Ice bag