

Snap25 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P60880

Item Number: CM0008201



Introduction

SNAP25 Rabbit Monoclonal Antibody (P60880) validated for WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, rat SNAP25, a t-SNARE protein involved in neurotransmitter release and synaptic function.

[Learn More](#)

Parameter	Value
Product Name	SNAP25 Rabbit mAb
Remarks/Alias	SUP; RIC4; SEC9; SNAP; CMS18; RIC-4; SNAP-25; bA416N4.2; dJ1068F16.2
Species	Human
Gene ID (Human)	6616
Gene ID	6616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human SNAP25 (NP_001309839.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P60880
Protein Weight	23 kDa
Shipping	Ice bag