

Snap91 Rabbit Polyclonal Antibody For Wb, Elisa - O60641

Item Number: CM0029399

Introduction

SNAP91 Rabbit pAb (O60641) is a rabbit polyclonal antibody recognizing clathrin coat assembly protein AP180 (CALM, AP180) in human, mouse, and rat. Validated for WB and ELISA, it aids vesicle trafficking research.

[Learn More](#)

Parameter	Value
Product Name	SNAP91 Rabbit pAb
Remarks/Alias	CALM, AP180, SNAP91
Species	Human
Gene ID	9892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 321-600 of human SNAP91 (NP_001229723.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60641
Protein Weight	93kDa
Shipping	Ice bag