

Snrpa1 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - P09661

Item Number: CM0002593



Introduction

SNRPA1 Rabbit Monoclonal Antibody (P09661) validated for WB, IP, IHC-P, ELISA in Human, Mouse, Rat. Core U2 snRNP component for pre-mRNA splicing. Ideal for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	SNRPA1 Rabbit mAb
Remarks/Alias	Lea1; U2A'; SNRPA1
Species	Human
GeneID	6627
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 156-255 of human SNRPA1 (P09661).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09661
Protein Weight	28 kDa
Shipping	Ice bag