

Snrpd2 Rabbit Pab For Wb, Ihc-P, Elisa - P62316

Item Number: CM0025197



Introduction

SNRPD2 Rabbit pAb is a polyclonal antibody targeting human SNRPD2 (Sm-D2), validated for WB, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for splicing research.

[Learn More](#)

Parameter	Value
Product Name	SNRPD2 Rabbit pAb
Remarks/Alias	SMD2; Sm-D2; SNRPD1; SNRPD2
Species	Human
GeneID	6633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human SNRPD2 (NP_004588.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62316
Protein Weight	14kDa
Shipping	Ice bag