

Snrpn Rabbit Polyclonal Antibody For Wb And Elisa - P63162

Item Number: CM0020215



Introduction

SNRPN rabbit polyclonal antibody validated for WB and ELISA applications. Detects human, mouse, and rat SNRPN, a splicing-associated nuclear protein involved in alternative RNA processing.

[Learn More](#)

Parameter	Value
Product Name	SNRPN Rabbit pAb
Remarks/Alias	SMN; PWCR; SM-D; sm-N; RT-LI; HCERN3; SNRNP-N; SNURF-SNRPN; SNRPN
Species	Human
Gene ID	6638
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SNRPN (NP_003088.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63162
Protein Weight	25kDa
Shipping	Ice bag