

Snx18 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Rf0

Item Number: CM0015922

Introduction

Rabbit polyclonal antibody targeting SNX18, validated for Western blot and ELISA applications. Reacts with human, mouse, and rat. Ideal for studying endocytosis, vesicle trafficking, and membrane dynamics. Specificity confirmed by recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	SNX18 Rabbit pAb
Remarks/Alias	SNAG1, SH3PX2, SH3PXD3B, SNX18
Species	Human
GeneID (Human)	112574
GeneID	112574
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 445-624 of human SNX18 (NP_001096045.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RF0
Protein Weight	69kDa
Shipping	Ice bag