

Snx4 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - 095219

Item Number: CM0003658

Introduction

High-quality SNX4 rabbit monoclonal antibody validated for Western blot, immunoprecipitation, and ELISA. Reacts with human and mouse SNX4, a sorting nexin involved in endocytosis and autophagy. Ideal for studying intracellular trafficking and receptor recycling.

[Learn More](#)

Parameter	Value
Product Name	SNX4 Rabbit mAb
Remarks/Alias	ATG24B; SNX4
Species	Human
Gene ID (Human)	8723
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95219
Protein Weight	52 kDa
Shipping	Ice bag