


Sod3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08294

Item Number: CM0026101



Introduction

SOD3 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat extracellular superoxide dismutase, a key antioxidant enzyme. Ideal for oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	SOD3 Rabbit pAb
Remarks/Alias	EC-SOD; SOD3
Species	Human
GeneID	6649
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 19-160 of human SOD3 (NP_003093.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08294
Protein Weight	26kDa
Shipping	Ice bag