

Somatostatin Receptor 2 (Sstr2) Rabbit Pab - P30874

Item Number: CM0025887



Introduction

The rabbit polyclonal antibody targeting Somatostatin Receptor 2 (SSTR2) is suitable for WB, IHC-P, ELISA; validated in human, mouse, and rat. It recognizes a GPCR involved in somatostatin-mediated hormone secretion inhibition.

[Learn More](#)

Parameter	Value
Product Name	Somatostatin Receptor 2 (SSTR2) Rabbit pAb
Remarks/Alias	SST2; Somatostatin Receptor 2 (SSTR2)
Species	Human
GeneID (Human)	6752
GeneID	6752
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-369 of human Somatostatin Receptor 2 (SSTR2) (NP_001041.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30874
Protein Weight	41kDa
Shipping	Ice bag