

Sox10 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P56693

Item Number: CM0007965

Introduction

Recombinant rabbit monoclonal antibody against human SOX10, validated for WB, IF/ICC and ELISA. Recognizes human, mouse, rat SOX10 (predicted 50 kDa). For neural crest, myelination, and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	SOX10 Rabbit mAb
Remarks/Alias	DOM; WS4; PCWH; WS2E; WS4C; SOX-10; SOX10
Species	Human
GeneID (Human)	6663
GeneID	6663
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56693
Protein Weight	50kDa
Shipping	Ice bag