

Sox10 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P56693

Item Number: CM0018970

Introduction

High-quality SOX10 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Targets transcription factor SOX-10, key for glial development and myelination. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	SOX10 Rabbit pAb
Remarks/Alias	DOM; WS4; PCWH; WS2E; WS4C; SOX-10; SOX10
Species	Human
Gene ID (Human)	6663
Gene ID	6663
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 270-360 of human SOX10 (NP_008872.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56693
Protein Weight	50kDa
Shipping	Ice bag