

Sox10 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P56693

Item Number: CM0030009

Introduction

Anti-SOX10 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat SOX10, a transcription factor crucial for oligodendrocyte maturation and myelination. Ideal for neuroscience and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	SOX10 Rabbit pAb
Remarks/Alias	DOM; WS4; PCWH; WS2E; WS4C; SOX-10; SOX10
Species	Human
GeneID (Human)	6663
GeneID	6663
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 270-360 of human SOX10 (NP_008872.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56693
Protein Weight	50kDa
Shipping	Ice bag