


Sox11 Rabbit Pab For Wb, If/Icc, Elisa - P35716

Item Number: CM0020279



Introduction

SOX11 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat SOX11. Targets transcription factor SOX-11 (P35716) involved in neurogenesis and organogenesis.

[Learn More](#)

Parameter	Value
Product Name	SOX11 Rabbit pAb
Remarks/Alias	CSS9; MRD27; IDDMOH; SOX11
Species	Human
GeneID (Human)	6664
GeneID	6664
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SOX11 (NP_003099.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35716
Protein Weight	47kDa
Shipping	Ice bag