


Sox2 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P48431

Item Number: CM0014090



Introduction

Anti-SOX2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Designed for stem cell and developmental biology research; reacts with human, mouse, and rat SOX2.

[Learn More](#)

Parameter	Value
Product Name	SOX2 Rabbit pAb
Remarks/Alias	ANOP3; MCOPS3; SOX2
Species	Human
GeneID (Human)	6657
GeneID	6657
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human SOX2 (NP_003097.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48431
Protein Weight	34kDa
Shipping	Ice bag