

Sox9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P48436

Item Number: CM0022002



Introduction

SOX9 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets transcription factor SOX9 essential for chondrogenesis and skeletal development. Suitable for IVD research.

[Learn More](#)

Parameter	Value
Product Name	SOX9 Rabbit pAb
Remarks/Alias	CMD1; SRA1; CMPD1; SRXX2; SRXY10; SOX9
Species	Human
Gene ID (Human)	6662
Gene ID	6662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-339 of human SOX9 (NP_000337.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48436
Protein Weight	56kDa
Shipping	Ice bag