

Sp140 Rabbit Pab For Wb And Elisa - Q13342

Item Number: CM0018790



Introduction

Rabbit polyclonal antibody against human SP140 (LYSp100), validated for Western blot and ELISA, cross-reacts with mouse and rat, detects 98kDa nuclear body protein.

[Learn More](#)

Parameter	Value
Product Name	SP140 Rabbit pAb
Remarks/Alias	LYSP100; LYSP100-A; LYSP100-B; SP140
Species	Human
GeneID (Human)	11262
GeneID	11262
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SP140 (NP_001265380)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13342
Protein Weight	98kDa
Shipping	Ice bag