

Sp4 Rabbit Polyclonal Antibody For Wb And Elisa - Q02446

Item Number: CM0024357

Introduction

Rabbit polyclonal antibody targeting human SP4 transcription factor, validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for transcriptional regulation and brain research.

[Learn More](#)

Parameter	Value
Product Name	SP4 Rabbit pAb
Remarks/Alias	HF1B; SPR-1; SP4
Species	Human
GeneID (Human)	6671
GeneID	6671
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-440 of human SP4 (NP_003103.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02446
Protein Weight	82kDa
Shipping	Ice bag