

Sparc Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09486

Item Number: CM0013153



Introduction

Rabbit polyclonal antibody targeting human SPARC (Osteonectin, BM-40). Validated for WB, IHC-P, IF/ICC, and ELISA applications with cross-reactivity to mouse and rat. Ideal for extracellular matrix and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SPARC Rabbit pAb
Remarks/Alias	ON; ONT; OI17; BM-40; SPARC
Species	Human
Gene ID (Human)	6678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-303 of human SPARC (NP_003109.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09486
Protein Weight	35kDa
Shipping	Ice bag