

Spred1 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z699

Item Number: CM0013625



Introduction

SPRED1 rabbit polyclonal antibody for WB and ELISA, cross-reactive with mouse. Targets human SPRED1, a tyrosine kinase substrate inhibiting MAP kinase activation, suitable for cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	SPRED1 Rabbit pAb
Remarks/Alias	LGSS; NFLS; hSpred1; spred-1; PPP1R147; SPRED1
Species	Human
GeneID (Human)	161742
GeneID	161742
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-320 of human SPRED1 (NP_689807.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q7Z699
Protein Weight	50kDa
Shipping	Ice bag