


Sptlc2 Rabbit Pab For Wb, Elisa - O15270

Item Number: CM0014591



Introduction

SPTLC2 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human, mouse, and rat serine palmitoyltransferase 2, a key enzyme in sphingolipid biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	SPTLC2 Rabbit pAb
Remarks/Alias	LCB2; SPT2; HSN1C; LCB2A; NSAN1C; hLCB2a; SPTLC2
Species	Human
Gene ID (Human)	9517
Gene ID	9517
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 423-562 of human SPTLC2 (NP_004854.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15270
Protein Weight	63kDa
Shipping	Ice bag