

Srbd1 Rabbit Pab - Q8N5C6

Item Number: CM0018149

Introduction

High-quality SRBD1 Rabbit Polyclonal Antibody validated for WB and ELISA, targeting human, mouse, and rat SRBD1 (UniProt Q8N5C6). Ideal for chromosome segregation, mitosis, and cell senescence research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | SRBD1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SRBD1                                                                                                              |
| Species          | Human                                                                                                              |
| GeneID (human)   | 55133                                                                                                              |
| GeneID           | 55133                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 330-600 of human SRBD1 (NP_060549.4) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8N5C6                                                                                                             |
| Protein Weight   | 112kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |