

Srek1 Rabbit Pab - Q8Wxa9

Item Number: CM0014238

Introduction

SREK1 Rabbit Polyclonal Antibody validated for WB and ELISA applications. Cross-reacts with human and mouse SREK1, a splicing regulatory protein involved in alternative splicing regulation. Suitable for spliceosome research.

[Learn More](#)

Parameter	Value
Product Name	SREK1 Rabbit pAb
Remarks/Alias	SFRS12; SRrp86; SRrp508; SREK1
Species	Human
Gene ID (Human)	140890
Gene ID	140890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SREK1 (NP_631907.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WXA9
Protein Weight	59kDa
Shipping	Ice bag