

Srsf10 Rabbit Pab For Wb And Elisa - 075494

Item Number: CM0026573

Introduction

Rabbit polyclonal antibody against human SRSF10 (Serine/arginine-rich splicing factor 10, SRp38). Validated for WB and ELISA; reacts with human, mouse, and rat. Ideal for studying pre-mRNA splicing, alternative exon skipping, and phosphorylation-regulated splicing repression.

[Learn More](#)

Parameter	Value
Product Name	SRSF10 Rabbit pAb
Remarks/Alias	NSSR, TASR, SRp38, TASR1, TASR2, FUSIP1, FUSIP2, SFRS13, SRp40, SFRS13A, PPP1R149, SRSF10
Species	Human
GeneID (Human)	10772
Immunogen	Synthetic peptide corresponding to amino acids 163-262 of human SRSF10 (NP_473357.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75494
Protein Weight	31kDa
Shipping	Ice bag