

Ssc4D Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wtu2

Item Number: CM0019140

Introduction

SSC4D Rabbit Polyclonal Antibody for WB and ELISA. Detects human, mouse, rat SSC4D, a secreted member of the scavenger receptor cysteine-rich superfamily. Suitable for studying its involvement in immune processes and proteomics.

[Learn More](#)

Specification	Detail
Product Name	SSC4D Rabbit pAb
Remarks/Alias	SRCRB4D, S4D-SRCRB, SRCRB-S4D, SSC4D
Species	Human
GeneID (Human)	136853
GeneID	136853
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 446-575 of human SSC4D (NP_542782.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTU2
Protein Weight	61kDa
Shipping	Ice bag