


Ssh2 Rabbit Polyclonal Antibody For Wb And Elisa - Q76I76

Item Number: CM0030831



Introduction

High-quality rabbit polyclonal antibody targeting human SSH2 protein. Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying actin dynamics and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	SSH2 Rabbit pAb
Remarks/Alias	SSH-2; SSH-2L; SSH2
Species	Human
GeneID (Human)	85464
GeneID	85464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 53-202 of human SSH2 (NP_203747.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q76I76
Protein Weight	158kDa
Shipping	Ice bag