

Ssx5 Rabbit Polyclonal Antibody For Wb, Elisa - O60225

Item Number: CM0024871

Introduction

Rabbit polyclonal antibody targeting human Protein SSX5 (O60225), validated for Western blot and ELISA. Cross-reacts with mouse and rat, detecting a 22 kDa band. Supports transcription regulation studies.

[Learn More](#)

Parameter	Description
Product Name	SSX5 Rabbit pAb
Remarks/Alias	SSX5
Species	Human
Gene ID (Human)	6758
Gene ID	6758
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-229 of human SSX5 (NP_066295.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60225
Protein Weight	22kDa
Shipping	Ice bag