

St3Gal3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q11203

Item Number: CM0025896

Introduction

Rabbit polyclonal antibody to ST3GAL3, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat. Useful for studying sialyltransferase function and glycosylation in Golgi.

[Learn More](#)

Parameter	Value
Product Name	ST3GAL3 Rabbit pAb
Remarks/Alias	ST3N; DEE15; MRT12; SIAT6; EIEE15; ST3GALII; ST3GalIII; ST3Gal III; ST3GAL3
Species	Human
Gene ID (Human)	6487
Gene ID	6487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-375 of human ST3GAL3 (NP_006270.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q11203
Protein Weight	42kDa
Shipping	Ice bag