

Stat1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P42224

Item Number: CM0013611



Introduction

Rabbit polyclonal antibody targeting human STAT1, validated for WB, IHC-P, IF/ICC, IP, ELISA, with cross-reactivity to mouse and rat. Ideal for interferon signaling, transcription regulation, and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	STAT1 Rabbit pAb
Remarks/Alias	CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91; STAT1
Species	Human
Gene ID (Human)	6772
Gene ID	6772
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 650-750 of human STAT1 (NP_009330.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42224
Protein Weight	87kDa
Shipping	Ice bag