

Stat3 Rabbit Pab For Wb, If/Icc, Elisa - P40763

Item Number: CM0015034



Introduction

High-quality STAT3 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat STAT3. Suitable for studying JAK/STAT signaling, inflammation, and cancer.

[Learn More](#)

Parameter	Value
Product Name	STAT3 Rabbit pAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1; STAT3
Species	Human
GeneID (Human)	6774
GeneID	6774
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human STAT3 (NP_003141.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40763
Protein Weight	88kDa
Shipping	Ice bag