

Stat3 Rabbit Pab For Wb, If/Icc, Ihc-P, Elisa - P40763

Item Number: CM0015936



Introduction

STAT3 Rabbit pAb (P40763) is a polyclonal antibody validated for WB, IF/ICC, IHC-P, ELISA in human, mouse, rat. It targets signal transducer and activator of transcription 3, ideal for cytokine signaling and immune regulation research.

[Learn More](#)

Parameter	Value
Product Name	STAT3 Rabbit pAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1; STAT3
Species	Human
GeneID (Human)	6774
GeneID	6774
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human STAT3 (NP_003141.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40763
Protein Weight	88 kDa
Shipping	Ice bag