


Stat5 Rabbit Monoclonal Antibody For Wb And Elisa - P42229/P51692

Item Number: CM0012527



Introduction

High-quality Anti-STAT5 Rabbit Monoclonal Antibody for Western blot and ELISA, recognizing human, mouse, and rat STAT5A/B. Ideal for studying cytokine signaling, transcription regulation, and lactation biology.

[Learn More](#)

Parameter	Value
Product Name	STAT5 Rabbit mAb
Remarks/Alias	MGF; STAT5
Species	Human
Gene ID (Human)	6776/6777
Gene ID	6776/6777
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-786 of human STAT5 (P42229)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42229/P51692
Protein Weight	90kDa
Shipping	Ice bag