


Stat5A Rabbit Mab For Wb, Elisa - P42229

Item Number: CM0002632



Introduction

STAT5A Rabbit Monoclonal Antibody validated for WB and ELISA. Recognizes human and mouse STAT5A, a transcription factor involved in cytokine signaling. Ideal for research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	STAT5A Rabbit mAb
Remarks/Alias	MGF; STAT5; STAT5A
Species	Human
GeneID (Human)	6776
GeneID	6776
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42229
Protein Weight	91kDa
Shipping	Ice bag