

Staufen Rabbit Monoclonal Antibody For Wb, Elisa - 095793

Item Number: CM0009883



Introduction

Staufen rabbit monoclonal antibody validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying RNA-binding proteins, mRNA localization, and viral host interactions. Full product specs available.

[Learn More](#)

Parameter	Value
Product Name	Staufen Rabbit mAb
Remarks/Alias	STAU; PPP1R150; Staufen
Species	Human
Gene ID (Human)	6780
Gene ID	6780
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 478-577 of human Staufen (095793)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095793
Protein Weight	63kDa
Shipping	Ice bag