


Steap3 Rabbit Polyclonal Antibody For Wb And Elisa - Q658P3

Item Number: CM0014350



Introduction

High-quality rabbit polyclonal antibody against human STEAP3, validated for WB and ELISA. Suitable for detecting metalloredutase STEAP3 in human and mouse samples, with a calculated molecular weight of 55 kDa.

[Learn More](#)

Parameter	Value
Product Name	STEAP3 Rabbit pAb
Remarks/Alias	STMP3; TSAP6; pHyde; AHMIO2; dudlin-2; dudulin-2; STEAP3
Species	Human
GeneID (Human)	55240
GeneID	55240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human STEAP3 (NP_060704.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q658P3
Protein Weight	55kDa
Shipping	Ice bag