

Stim1 Rabbit Monoclonal Antibody For Wb, Elisa - Q13586

Item Number: CM0010626



Introduction

STIM1 Rabbit mAb is a recombinant monoclonal antibody validated for Western blot and ELISA, reacting with human, mouse, and rat STIM1. Ideal for studying calcium signaling and store-operated calcium entry.

[Learn More](#)

Parameter	Value
Product Name	STIM1 Rabbit mAb
Remarks/Alias	GOK; TAM; TAM1; IMD10; STRMK; D11S4896E; STIM1
Species	Human
Gene ID (Human)	6786
Gene ID	6786
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 586-685 of human STIM1 (Q13586).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13586
Protein Weight	77kDa
Shipping	Ice bag