

Stip1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P31948

Item Number: CM0017822



Introduction

Rabbit polyclonal antibody against STIP1 (HOP), validated for WB, IHC-P, ELISA in human, mouse, rat. Co-chaperone mediating HSP70-HSP90 interaction. Ideal for stress response and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	STIP1 Rabbit pAb
Remarks/Alias	HOP, P60, STI1, STI1L, HEL-S-94n, IEF-SSP-3521, STIP1
Species	Human
GeneID (Human)	10963
GeneID	10963
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human STIP1 (NP_006810.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31948
Protein Weight	63kDa
Shipping	Ice bag