

Stoml2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ujz1

Item Number: CM0017685

Introduction

STOML2 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat STOML2 (39kDa, Q9UJZ1). Ideal for mitochondrial research, T-cell activation, and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	STOML2 Rabbit pAb
Remarks/Alias	SLP-2; HSPC108; STOML2
Species	Human
GeneID (Human)	30968
GeneID	30968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-356 of human STOML2 (NP_038470.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJZ1
Protein Weight	39kDa
Shipping	Ice bag